


MAX Dental Implants



*Note: No flute in last two threads to protect sinus

* The Tri-MAX implant of 7.0 mm lengths features a 1.6 mm smooth collar to accommodate the internal connection.



Internal Tri-Lobe



Internal Octagon/Morse Taper



External Hexagon

Tri-MAX:		
DIAMETER	LENGTH	PLATFORM
7.0 mm	7.0 mm	Ø5.0 mm
7.0 mm	9.0 mm	Ø5.0 mm
7.0 mm	11.0 mm	Ø5.0 mm
8.0 mm	7.0 mm	Ø6.0 mm
8.0 mm	9.0 mm	Ø6.0 mm
8.0 mm	11.0 mm	Ø6.0 mm
9.0 mm	7.0 mm	Ø6.0 mm
9.0 mm	9.0 mm	Ø6.0 mm
9.0 mm	11.0 mm	Ø6.0 mm

MAX-IT:		
DIAMETER	LENGTH	PLATFORM
7.0 mm	7.0 mm	Ø4.8 mm
7.0 mm	9.0 mm	Ø4.8 mm
7.0 mm	11.0 mm	Ø4.8 mm
8.0 mm	7.0 mm	Ø6.5 mm
8.0 mm	9.0 mm	Ø6.5 mm
8.0 mm	11.0 mm	Ø6.5 mm
9.0 mm	7.0 mm	Ø6.5 mm
9.0 mm	9.0 mm	Ø6.5 mm
9.0 mm	11.0 mm	Ø6.5 mm

MAX:		
DIAMETER	LENGTH	PLATFORM
7.0 mm	7.0 mm	Ø5.0 mm
7.0 mm	9.0 mm	Ø5.0 mm
7.0 mm	11.0 mm	Ø5.0 mm
8.0 mm	7.0 mm	Ø6.0 mm
8.0 mm	9.0 mm	Ø6.0 mm
8.0 mm	11.0 mm	Ø6.0 mm
9.0 mm	7.0 mm	Ø7.0 mm
9.0 mm	9.0 mm	Ø7.0 mm
9.0 mm	11.0 mm	Ø7.0 mm

SURGICAL AND PROSTHETIC KEY BENEFITS

- Immediate placement in a molar extraction socket
- Maintains the natural architecture of the site and maximizes bone preservation
- Restoration 4-months or less after placement
- Proven thread pattern allows for superior primary stability
- Industry compatible connection options
- Allows for platform switching



5 Holland Building 209 Irvine, CA 92618
Toll Free 866 700-2100 **Tel** 949 273-8505 **Fax** 949 273-8508
www.southernimplants.us
 ©2011 Southern Implants Inc. All rights reserved. SIM-11-16MX



Bright Ideas in Dental Implants
Since 1987